

DEVICE PRESSURE SETTINGS

THE BREATHER.

THE BREATHER - RESPIRATORY IMPAIRED ADULTS			
Full Range (impaired, 30 lpm): -52 to +31 cm H ₂ O			
INHALE		EXHALE	
SETTING	PRESSURE	SETTING	PRESSURE
1	-10 cm H ₂ O	1	5 cm H ₂ O
2	-12 cm H ₂ O	2	13 cm H ₂ O
3	-15 cm H ₂ O	3	21 cm H ₂ O
4	-28 cm H ₂ O	4	23 cm H ₂ O
5	-38 cm H ₂ O	5	31 cm H ₂ O
6	-52 cm H ₂ O		
Flow Rate: 30 lpm		Flow Rate: 30 lpm	

THE BREATHER - HEALTHY ADULTS			
Full Range (healthy, 70 lpm): -243 to +126 cm H ₂ O			
INHALE		EXHALE	
SETTING	PRESSURE	SETTING	PRESSURE
1	-35 cm H ₂ O	1	20 cm H ₂ O
2	-46 cm H ₂ O	2	53 cm H ₂ O
3	-64 cm H ₂ O	3	80 cm H ₂ O
4	-117 cm H ₂ O	4	107 cm H ₂ O
5	-180 cm H ₂ O	5	126 cm H ₂ O
6	-243 cm H ₂ O		
Flow Rate: 60 lpm		Flow Rate: 70 lpm	

BREATHER. FIT

BREATHER FIT - HEALTHY ADULTS & ATHLETES			
Full Range (healthy, 70 lpm): -335 to +250 cm H ₂ O			
INHALE		EXHALE	
SETTING	PRESSURE	SETTING	PRESSURE
1	-64 cm H ₂ O	1	29 cm H ₂ O
2	-74 cm H ₂ O	2	71 cm H ₂ O
3	-89 cm H ₂ O	3	93 cm H ₂ O
4	-145 cm H ₂ O	4	135 cm H ₂ O
5	-240 cm H ₂ O	5	250 cm H ₂ O
6	-335 cm H ₂ O		
Flow Rate: 60 lpm		Flow Rate: 70 lpm	

Additional technical specifications of The Breather www.pnmedical.com/product/the-breather/

Additional technical specifications of Breather Fit www.pnmedical.com/product/breather-fit/

CONSIDERATIONS

Indicated low and pressure ranges may vary according to respiratory capacity and should be considered as guidelines.